



ELK Biotechnology

 15000余种抗体

 WB、IHC、IF

 常用多抗单抗二抗

 一百种内参抗体

 文献引用数量超 3000篇

集研发、生产、营销一体化全球化试剂供应商!!!

高分文献

27.4

Advanced materials

优选抗体



Disrupting Intracellular Homeostasis by Copper-Based Nanoinducer with Multiple Enzyme-Mimicking Activities to Induce Disulfidptosis-Enhanced Pyroptosis for Tumor Immunotherapy

DOI: <https://doi.org/10.1002/adma.202410957>

引用产品&名称: ES7676 Cleaved-Caspase-1 (M211) rabbit pAb

公司简介

科鹿生物科技公司位于武汉精准医疗产业基地。ELK生物科技公司成立于2018年，是一家专业从事生命科学和免疫研究的高科技生物企业。公司拥有3000平方米的实验室，拥有独立的细胞培养室和SPF动物室，专业从事形态病理学、细胞生物学、分子生物学、蛋白质免疫学、模型动物相关产品的引进和研发。公司成立以来，团队成员参与了多项医学生物学研究项目和国家自然科学基金项目。并与全国多家企业和科研机构建立合作关系，为每一位科研人员提供最优质、最高效的服务和产品。



科鹿（武汉）生物科技有限责任公司

☎ +86-027-59760950 ✉ elkbio@elkbiotech.com

📍 武汉东湖新技术开发区高科园三路9号光谷精准医疗产业基地2.1期11栋



为什么选择ELK抗体？

- 1 多重验证、特异性高、反应物种广、稀释比例佳
- 2 一致性高、可重复性高，高灵敏度、超强特异性
- 3 验证实验类型ELISA、WB、IHC、IF
- 4 严格的质控流程 100%无风险保障、无忧售后
- 5 价格一致50ul/1280元，HRP二抗100ul/160元




抗体
推荐

ELK BIOTECHNOLOGY


ELK 可提供经过WB, IHC, IF等实验验证的高品质抗体 15000余种，覆盖人类基因组靶标。涉及生命科学基础代谢与转化研究的主要领域，包括表观遗传学和转录调节、代谢与内分泌、信号转导、细胞生物学、心血管、神经科学、干细胞和发育、肿瘤抑制和细胞凋亡、炎症与免疫、蛋白质翻译后修饰等。

产品货号	产品名称	产品货号	产品名称
EA152-2	E-Cadherin Rabbit pAb	EM1207-2	STAT3 Mouse mAb
SA003-1	HRP-conjugated Goat Anti-Rabbit IgG	SA002-2	HRP-conjugated Goat Anti-Mouse IgG
EA169-2	GSK 3β Rabbit pAb	ES1255-2	Akt (phospho Ser473) Rabbit pAb
EM1159-2	Vimentin Mouse mAb	EM1111-0	NFκB p65 Mouse mAb
EM1123-2	Collagen III Mouse mAb	ES1530-2	ERK 1/2 (phospho Thr202/Y204) Rabbit pAb
SA002-1	HRP-conjugated Goat Anti-Mouse IgG	EA014-1	β-actin Rabbit pAb
EM1022-1	His-Tag(HRP Conjugated) Mouse mAb	ES1442-2	NOS3 (phospho Ser1177) Rabbit pAb
EA234-2	STAT3 Rabbit pAb	SA003-2	HRP-conjugated Goat Anti-Rabbit IgG
EM1137-2	Fibronectin Mouse mAb	EM1066-2	NSE Mouse mAb




EM1052

P53单克隆抗体 HeLa细胞免疫荧光染色




EM1029


β-Tubulin 单克隆抗体 HeLa细胞免疫荧光染色



GAPDH Western blot analysis of HeLa (1), Rat brain (2), Rabbit Muscle (3), Sheep Muscle (4), and Mouse brain (5) with GAPDH mouse mAb (2B8) diluted at 110,000.




Western blot analysis of HepG2 (1), Rat liver (2), Mouse kidney (3), Rabbit testis (4), Sheep lung (5), 293T (6) with β-Actin rabbit polyclonal antibody diluted at 15000.



EA337

HSP40 多克隆抗体
结肠癌组织免疫组化染色



EM1061

CEA单克隆抗体
结肠癌组织免疫组化染色