

Collagen III Mouse mAb Catalog NO.: EM1123 For research use only.

## Overview

Product name Collagen III Mouse Monoclonal antibody

**Source** Mouse

Applications WB IHC

Species reactivity Human Rat Mouse

Recommended dilutions WesternBlot:1/1000

Immunohistochemistry:1/200-500

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

**Clonality** Monoclonal

Concentration 1 mg/ml

Observed band 138kDa

GenelD (Human) 1281

Human Swiss-Prot No. P02461

**Cellular localization** Secretedextracellular spaceextracellular matrix


Alternative Names CO3A1 antibody Collagen alpha(III) chain antibody Collagen III alpha chain

precursor antibody Collagen type III alpha antibody

**Background** Collagen alpha-1(III) chain is a fibrillar collagen that is found in extensible

connective tissues such as skin lung and the vascular system frequently in

association with type I collagen.



Western blot analysis of Hela with EM1123 diluted at:1000.

Immunohistochemical analysis of paraffin-embedded human Skin Tissue using Collagen III (EM1123) Mouse mAb diluted at:500.