

p53 (Acetyl Lys372) rabbit pAb

Cat No.:ES1111

For research use only

Overview

Immunogen

Product Name p53 (Acetyl Lys372) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300

ELISA: 1/20000. Not yet tested in other applications. The antiserum was produced against synthesized

Acetyl-peptide derived from human p53 around the

Acetylation site of Lys372. AA range:331-380

Specificity Acetyl-p53 (K372) Polyclonal Antibody detects

endogenous levels of p53 protein only when

acetylated at K372.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Cellular tumor antigen p53

Gene Name TP53

Cellular localization Cytoplasm . Nucleus . Nucleus, PML body .

Endoplasmic reticulum . Mitochondrion matrix . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Recruited into PML bodies together with CHEK2 (PubMed:12810724).

Translocates to mitochondria upon oxidative stress (PubMed:22726440). Translocates to mitochondria

in response to mitomycin C treatment

(PubMed:27323408). .; [Isoform 1]: Nucleus . Cytoplasm. Predominantly nuclear but localizes to the cytoplasm when expressed with isoform 4.; [Isoform 2]: Nucleus. Cytoplasm. Localized mainly in the nucleus with minor staining in the cytoplasm.; [Isoform 3]: Nucleus. Cytoplasm. Localized in the nucleus in most cells but found in the cytoplasm in

some cells.; [Isoform 4]: Nucleus. Cytoplasm.

+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com





Purification

Predominantly nuclear but translocates to the cy The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 53kD
Human Gene ID 7157
Human Swiss-Prot Number P04637


Alternative Names TP53:

TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor

p53

Background


tumor protein p53(TP53) Homo sapiens This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons (PMIDs: 12032546, 20937277). [provided by RefSeq, Feb 2013],




Western Blot analysis of Jurkat cells using Acetyl-p53 (K372) Polyclonal Antibody. Antibody was diluted at 1:2000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000








Western Blot analysis of 293T cells using Acetyl-p53 (K372) Polyclonal Antibody diluted at 1:2000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

Immunohistochemical analysis of paraffin-embedded human-breast, antibody was diluted at 1:100



Western blot analysis of lysate from Jurkat cells, using p53 (Acetyl-Lys372) Antibody.

