

## Bcl-2 (phospho Ser70) rabbit pAb

## Cat No.:ES1272

For research use only

## Overview

Product Name	Bcl-2 (phospho Ser70) rabbit pAb	
Host species	Rabbit	
Applications	WB;IHC;IF;IP;ELISA	
Species Cross-Reactivity	Human;Monkey	
Recommended dilutions	Western Blot: 1/500 - 1/2000.	
	Immunohistochemistry: 1/100 - 1/300.	
	Immunoprecipitation: 2-5 ug/mg lysate. ELISA:	
	1/10000. Not yet tested in other applications.	
Immunogen	The antiserum was produced against synthesized	
	peptide derived from human BCL-2 around the	
	phosphorylation site of Ser70. AA range:41-90	
Specificity	Phospho-Bcl-2 (S70) Polyclonal Antibody detects	
	endogenous levels of Bcl-2 protein only when	
	phosphorylated at S70.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Apoptosis regulator Bcl-2	
Gene Name	BCL2	
<b>Cellular localization</b>	Mitochondrion outer membrane ; Single-pass	
	membrane protein . Nucleus membrane ;	
	Single-pass membrane protein . Endoplasmic	
	reticulum membrane ; Single-pass membrane	
	protein . Cytoplasm .	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	llı.
Observed band	26kD	
Human Gene ID	596	
Human Swiss-Prot Number	P10415	
Alternative Names	BCL2; Apoptosis regulator Bcl-2	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

VEC

Phospho-Bcl-2 (S70)

HeLa

(kD)

117-85-

48-

34-

26-

19-

158

100-70-55-40-35--

25-

15

BCL2, apoptosis regulator(BCL2) Homo sapiens This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],

Western Blot analysis of various cells using Phospho-Bcl-2 (S70) Polyclonal Antibody diluted at 1:1000

Western Blot analysis of HeLa cells using Phospho-Bcl-2 (S70) Polyclonal Antibody diluted at 1:1000

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using BCL-2 (Phospho-Ser70) Antibody. The picture on the right is blocked with the phospho peptide.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Western blot analysis of lysates from HeLa cells treated with LPS (40nM, 30mins), K562 cells treated with calyculin A (50ng/ml, 30mins) and COS-7 cells treated with H2O2 (1ng/ml, 15mins), using BCL-2 (Phospho-Ser70) Antibody. The lane on the right is blocked with the phospho peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C