

## Cdc25B (phospho Ser323) rabbit pAb

## Cat No.:ES1283

For research use only

## Overview

Product Name	Cdc25B (phospho Ser323) rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human CDC25B around the
	phosphorylation site of Ser323. AA range:289-338
Specificity	Phospho-Cdc25B (S323) Polyclonal Antibody detects
	endogenous levels of Cdc25B protein only when
	phosphorylated at S323.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	M-phase inducer phosphatase 2
Gene Name	CDC25B
<b>Cellular localization</b>	Cytoplasm, cytoskeleton, microtubule organizing
	center, centrosome . Cytoplasm, cytoskeleton,
	spindle pole .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70kD
Human Gene ID	994
Human Swiss-Prot Number	P30305
Alternative Names	CDC25B; CDC25HU2; M-phase inducer phosphatase
	2; Dual specificity phosphatase Cdc25B
Background	cell division cycle 25B(CDC25B) Homo sapiens
	CDC25B is a member of the CDC25 family of
	phosphatases. CDC25B activates the cyclin



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



dependent kinase CDC2 by removing two phosphate groups and it is required for entry into mitosis. CDC25B shuttles between the nucleus and the cytoplasm due to nuclear localization and nuclear export signals. The protein is nuclear in the M and G1 phases of the cell cycle and moves to the cytoplasm during S and G2. CDC25B has oncogenic properties, although its role in tumor formation has not been determined. Multiple transcript variants for this gene exist. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells usingPhospho-Cdc25B (S323) Polyclonal Antibody diluted atPhospho-Cdc25B (S323) 1:1000

Western Blot analysis of VEC cells using Phospho-Cdc25B (S323) Polyclonal Antibody diluted at 1:1000



VEC 293 A549

138:

70--55-40-35--

15---

138---

70---55---40---35---25---

15----

VEC

p-Cdc25B (\$323)

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absor



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using CDC25B (Phospho-Ser323) Antibody



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C