

EphA2 (phospho Tyr588) rabbit pAb

Cat No.:ES1306

For research use only

Overview

Product Name	EphA2 (phospho Tyr588) rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not
	yet tested in other applications.
Immunogen	Synthesized phospho-peptide around the
	phosphorylation site of human EphA2 (phospho
	Tyr588)
Specificity	Phospho-EphA2 (Y588) Polyclonal Antibody detects
	endogenous levels of EphA2 protein only when
	phosphorylated at Y588.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ephrin type-A receptor 2
Gene Name	EPHA2
Cellular localization	Cell membrane ; Single-pass type I membrane
	protein . Cell projection, ruffle membrane ;
	Single-pass type I membrane protein . Cell
	projection, lamellipodium membrane ; Single-pass
	type I membrane protein . Cell junction, focal
	adhesion . Present at regio
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	130kD
Human Gene ID	1969
Human Swiss-Prot Number	P29317
Alternative Names	EPHA2; ECK; Ephrin type-A receptor 2; Epithelial cell
	kinase; Tyrosine-protein kinase receptor ECK



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

This gene belongs to the ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events, particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. This gene encodes a protein that binds ephrin-A ligands. Mutations in this gene are the cause of certain genetically-related cataract disorders.[provided by RefSeq, May 2010],





0.8

0.7

0.6

0.5

0.4 0.3 0.2 0.1 Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using EphA2 (Phospho-Tyr588) Antibody



+86-27-59760950

Phosphopeptide

Non-Phosphopeptide

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Western blot analysis of lysates from HepG2 cell, using phospho-EphA2 (Phospho-Tyr588) antibody.

-	95
-	72
_	55

-170



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C