

MEK-2 (phospho Thr394) rabbit pAb

Cat No.:ES1356

For research use only

Overview

Product Name	MEK-2 (phospho Thr394) rabbit pAb	
Host species	Rabbit	
Applications	WB;IHC;IF;IP;ELISA	
Species Cross-Reactivity	Human;Mouse;Rat	
Recommended dilutions	Western Blot: 1/500 - 1/2000.	
	Immunohistochemistry: 1/100 - 1/300.	
	Immunoprecipitation: 2-5 ug/mg lysate. ELISA:	
	1/10000. Not yet tested in other applications.	
Immunogen	The antiserum was produced against synthesized	
	peptide derived from human MEK2 around the	
	phosphorylation site of Thr394. AA range:261-310	
Specificity	Phospho-MEK-2 (T394) Polyclonal Antibody detects	
	endogenous levels of MEK-2 protein only when	
	phosphorylated at T394.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Dual specificity mitogen-activated protein kinase	
	kinase 2	
Gene Name	MAP2K2	
Cellular localization	Cytoplasm . Membrane ; Peripheral membrane	
	protein . Membrane localization is probably	
	regulated by its interaction with KSR1	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration Observed band	1 mg/ml	
Human Gene ID	5606	
Human Swiss-Prot Number	5606 P36507	
Alternative Names		
Alternative indiffes	MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP	
	milogen-activated protein kinase kinase 2, MAP	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

VEC 293

138-130-

70-55-

40-35-25-

15-

VEC

p-MEK-2 (T394)

138-70--55-

40---35----25---

15----

kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2

The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase is known to play a critical role in mitogen growth factor signal transduction. It phosphorylates and thus activates MAPK1/ERK2 and MAPK2/ERK3. The activation of this kinase itself is dependent on the Ser/Thr phosphorylation by MAP kinase kinase kinases. Mutations in this gene cause cardiofaciocutaneous syndrome (CFC syndrome), a disease characterized by heart defects, mental retardation, and distinctive facial features similar to those found in Noonan syndrome. The inhibition or degradation of this kinase is also found to be involved in the pathogenesis of Yersinia and anthrax. A pseudogene, which is located on chromosome 7, has been identified for this gene. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using Phospho-MEK-2 (T394) Polyclonal Antibody diluted at 1:2000

Western Blot analysis of VEC cells using Phospho-MEK-2 (T394) Polyclonal Antibody diluted at 1:2000

www.elkbiotech.com



+86-27-59760950 ELKbio@ELKbiotech.com 23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C

Phospho-MEK-2 (T394)



human breast carcinoma, using MEK2 (Phospho-Thr394) Antibody. The picture on the right is blocked with the phospho peptide.



Western blot analysis of lysates from ovary cancer, using MEK2 (Phospho-Thr394) Antibody. The lane on the right is blocked with the phospho peptide.

Immunohistochemistry analysis of paraffin-embedded



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C