

FITC-conjugated Goat Anti-Rabbit IgG

Catalog NO.: SA006

For research use only.

Overview

Product name FITC-conjugated Goat Anti-Rabbit IgG

Recommended dilutions 1:20-1:100 for most applications

NOTE: Optimal dilutions should be determined by the end user.

Size $100\mu L$

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20°C for one year. Avoid repeated freeze-

thaw cycles.

Note This product is for research use only,not for diagnostic or

therapeutic use.

https://www.elkbiotech.com T

T:86-27-59760950

E: elkbio@elkbiotech.com

Fluorophore FITC

Absorption Peak 492nm

Emission Peak 520nm